

User Story: As a clinician, I want to see other biometric data, so that I can better treat my patients. (5 points)

User Story: As a clinician, I want to see patient data appear on the dashboard in real time, so I can respond to my patients' reactions quickly. (8 points).


Real-time graphs of patient data on ClinicDashBoard page

We implemented real-time graphs of 10 data streams on ClinicDashboard page, Clinician can monitor patient data in detail.



User Story: As a clinician, I want to export patient data, so that other people can use the data. (5 points)

Press the Export Data button to collect patient data to the backend



Patient data received by backend

| id | date_time_recorded | patientid | sessionid | value |
|------|----------------------------|-----------|---------------------|--------|
| 26 | 2024-04-27 12:31:01.379000 | 1 | 1714221051236076000 | 0.8187 |
| 27 | 2024-04-27 12:31:01.628000 | 1 | 1714221051236076000 | 0.7999 |
| 28 | 2024-04-27 12:31:01.703000 | 1 | 1714221051236076000 | 0.7815 |
| 29 | 2024-04-27 12:31:01.805000 | 1 | 1714221051236076000 | 0.7724 |
| 30 | 2024-04-27 12:31:01.903000 | 1 | 1714221051236076000 | 0.7547 |
| 31 | 2024-04-27 12:31:02.030000 | 1 | 1714221051236076000 | 0.7459 |
| 32 | 2024-04-27 12:31:02.113000 | 1 | 1714221051236076000 | 0.7373 |
| 33 | 2024-04-27 12:31:02.208000 | 1 | 1714221051236076000 | 0.7288 |
| 34 | 2024-04-27 12:31:02.295000 | 1 | 1714221051236076000 | 0.7203 |
| NULL | NULL | NULL | NULL | NULL |

User Story: As a clinician, I want to be alerted when a patient's vitals exceed threshold, so that I do not overexpose my patients. (5 points)

Heart rate alert system

If the heart rate exceeds 100 beats per minute, the line turns red. Conversely, if the heart rate drops below 50 beats per minute, the line turns blue.

In addition to the line color, we have implemented a blinking effect to further emphasize abnormal heart rate values:


When the heart rate exceeds 100 beats per minute, the heart rate display blinks in red. Similarly, if the heart rate falls below 50 beats per minute, the display blinks in blue.

We also made an abnormal heart rate recorder, when abnormal heart rate appeared, the recorder will record the abnormal heart rate and the Time of occurrence, and show the message through message box

High heart rate alert:



Low Heart rate:



Abnormal Heart rate recorder:



User Story: As a clinician, I want to have the sensor automatically connect to Enhancite without any third-party programs, so that the program is easy to use. (8 points)

To automatically connect to Enhancite, we made a `udpport.py` python program to receive data automatically sent by Emotibit device from UDP port, and we made a `UdpPortControl.py` python program to control `udpport.py`. But Users still need to use the EmotiBit Oscilloscope to send data to the UDP port.

udpport.py:

```

udpport.py > udp_to_websocket
1  import socket
2  import asyncio
3  import websockets
4  websocket_connections = set()
5  async def udp_to_websocket():
6      udp_sock = socket.socket(socket.AF_INET, socket.SOCK_DGRAM)
7      udp_sock.bind(('0.0.0.0', 12346))
8      udp_sock.setblocking(False) # Make the socket non-blocking
9
10     # Start the WebSocket server
11     async with websockets.serve(handler, "localhost", 8765,
12 ping_timeout=30, close_timeout=10):
13         print("WebSocket server started...")
14         while True:
15             await asyncio.sleep(0.016) # Small delay to yield control to the loop
16             try:
17                 data, addr = udp_sock.recvfrom(2048)
18                 print(f"Received UDP data: {data} from {addr}")
19                 await forward_data(data)
20             except BlockingIOError:
21                 continue # No data available, continue the loop
22             except Exception as e:
23                 print(f"Unexpected error: {e}") # Log unexpected errors
24
25     async def handler(websocket, path):
26         websocket_connections.add(websocket)
27         try:
28             await websocket.wait_closed()
29         finally:
30             websocket_connections.remove(websocket)
31
32     async def forward_data(data):
33         # Convert data to string
34         data_str = data.decode()
35         # Create a list of websockets that failed to receive the data
36         failed_websockets = set()
37         # Create a copy of the set for iteration to prevent modification errors
38         connections = set(websocket_connections)
39         for websocket in connections:
40             if websocket.open:
41                 try:
42                     await websocket.send(data_str)
43                 except websockets.exceptions.ConnectionClosed:
44                     failed_websockets.add(websocket)
45         # Remove failed websockets from the original set
46         websocket_connections.difference_update(failed_websockets)
47
48 if __name__ == "__main__":
49     asyncio.run(udp_to_websocket())

```

UdpPortControl.py:

```

UdpPortControl.py > ...
1  from flask import Flask, jsonify
2  import subprocess
3  import os
4  import signal
5  import platform
6
7  # from flask_cors import CORS
8
9  app = Flask(__name__)
10
11 process = None
12
13 @app.route('/start-udp', methods=['POST'])
14 def start_udp():
15     global process
16     if process is None:
17         if platform.system() == 'Windows':
18             # you can input your own python.exe path to replace python, make sure the path is correct
19             process = subprocess.Popen(['c:/Users/jayye999/Desktop/PP1/Enhancite/.venv/Scripts/python.exe', './udpport.py'], creationflags=subprocess.CREATE_NEW_PROCESS_GROUP)
20         else:
21             # you can input your own python.exe path to replace python, make sure the path is correct
22             process = subprocess.Popen(['c:/Users/jayye999/Desktop/PP1/Enhancite/.venv/Scripts/python.exe', './udpport.py'], preexec_fn=os.setsid)
23         return jsonify({'message': 'UDP script started'}), 200
24     return jsonify({'message': 'UDP script already running'}), 400
25
26 @app.route('/stop-udp', methods=['POST'])
27 def stop_udp():
28     global process
29     if process:
30         if platform.system() == 'Windows':
31             os.kill(process.pid, signal.CTRL_BREAK_EVENT)
32         else:
33             os.killpg(os.getpgid(process.pid), signal.SIGTERM)
34         process = None
35     return jsonify({'message': 'UDP script stopped'}), 200
36     return jsonify({'message': 'No UDP script running'}), 400
37
38
39 if __name__ == '__main__':
40     app.run(host='0.0.0.0', port=5000)
41

```